

initial

CERAMIC SYSTEM



INITIAL COLOR COMBINATION CHART

INITIAL COLOR COMBINATION CHART																
Vita Classic Shade	A1	A2	A3	A3.5	A4	B1	B2	B3	B4	C1	C2	C3	C4	D2	D3	D4
Powder Opaque	A1	A2	A3	A3.5	A4	B1	B2	B3	B4	C1	C2	C3	C4	D2	D3	D4
Opacious Dentine	A1	A2	A3	A3.5	A4	B1	B2	B3	B4	C1	C2	C3	C4	D2	D3	D4
Fluorescence Dentin	91	93	93	93	93	91	92	92	92	91	93 + IN 51	93 + IN 51	93 + IN 51	91	93 + IN 51	93 + IN 51
Enamel	E58	E58	E59	E59	E60	E57	E59	E59	E59	E60	E59	E59	E60	E60	E59	E59
Inside (high chromatic dentine colors)	IN-41	IN-44	IN-44	IN-44	IN-42	IN-44 minimum	IN-44	IN-47	IN-46	IN-51 minimum	IN-51	IN-51	IN-51	IN-44	IN-44	IN-44 IN-51
Inside (incisal colors & occlusal colors)	IN-41	IN-44	IN-44	IN-44	IN-42	IN-44 minimum	IN-44	IN-47	IN-46	IN-51 minimum	IN-51	IN-51	IN-51	IN-44	IN-44	IN-44 IN-51
Inside (cervical)	IN-41	IN-44	IN-44	IN-44	IN-42	IN-44 minimum	IN-44	IN-47	IN-46	IN-51 minimum	IN-51	IN-51	IN-51	IN-44	IN-44	IN-44 IN-51
Shoulder Opaque																
Shoulder Opaque	SO-37	SO-37	SO-37	SO-37	SO-37	SO-37	SO-37	SO-37	SO-37	SO-37	SO-37	SO-37	SO-37	SO-37	SO-37	SO-37
Shoulder Transparent	ST-31	ST-31	ST-31	ST-31	ST-36	ST-30	ST-30	ST-32	ST-32	ST-30	ST-30	ST-30	ST-35	ST-31	ST-31	
		90%	80%	50%		50%	20%		50%	50%	30%	20%		80%	50%	
		ST-36	ST-36	ST-36		ST-32	ST-32		ST-33	ST-35	ST-35	ST-35		ST-35	ST-35	
		10%	20%	50%		50%	80%		50%	50%	70%	81%		30%	50%	